
wellion®

SENSOR

DE Systembeschreibung

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System Description

Thank you for choosing the Wellion SENSOR Continuous Glucose Monitoring System (rtCGMS) (hereafter referred to as CGM). The Wellion SENSOR CGM System includes products such as Portable Diabetes Assistant, a Transmitter and a Glucose Sensor.

The Wellion SENSOR CGM provides real-time glucose levels and allows you to continuously view your sensor glucose values on your selected mobile device. The system tracks your glucose every minute by measuring the amount of glucose in the interstitial fluid. A sensor inserted in your skin, sends glucose results to the Portable Diabetes Assistant.

Note: Please read all the instructions provided in this Instruction for Use before using these products.

EN

Content

System Description 129
Important safety information 132
Indications for use 132
Intended Purpose 132
Indications 133
Patients 133
Intended User 133
Contraindications 134
Warning 134
Precautions 137
Potential clinical side-effects 138
Additional security information 138
Product components 139
Component description 142
Portable Diabetes Assistant 142
Wellion SENSOR Transmitter 142
Wellion SENSOR CGM package 142
Using a new glucose sensor 143
Applying a sensor 143
Pairing with a new transmitter 146
Sensor start-up 146

Sensor calibration 147
Removing a sensor 148
Replacing the sensor 149
Wellion AiDEX SENSOR App – Diabetes Management System 149
Caring for your Continuous Glucose Monitoring System (rtCGMS) 149
Cleaning the transmitter 149
Disposal 150
Transportation 150
Storage 150
Specifications 151
Electromagnetic compatibility 152
Clinical performance 157
Appendix 162
Symbols 162
Potential interference information 164
Potential Risks 164
Inaccurate glucose values 164
Mild to severe to sensor related-wear reactions 164
Hyperglycemia or hypoglycemia 165
Potential clinical benefits 165
Glossary 166

Important safety information

Indications for use

EN

The Continuous Glucose Monitoring System is a real- time, continuous glucose monitoring device. When the system is used together with compatible devices, it is indicated for the management of diabetes in people age 2 and older. It is designed to replace finger stick blood glucose testing for diabetes treatment decisions.

Interpretation of the system results should be based on the glucose trends and several sequential readings over time. The system also detects trends and tracks patterns, and aids in the detection of episodes of hyperglycemia and hypoglycemia, facilitating both acute and long-term therapy adjustment.

Intended Purpose

Continuous Glucose Monitoring System Sensor: When the Continuous Glucose Monitoring System Sensor is used together with compatible devices, it is intended to continuously measure the glucose in the interstitial fluid and is designed to replace fingerstick blood glucose (BG) testing for treatment decisions.

Continuous Glucose Monitoring System Transmitter: When the Continuous Glucose Monitoring System Transmitter is used together with compatible devices, it is intended to continuously measure the

glucose in the interstitial fluid and is designed to replace fingerstick blood glucose (BG) testing for treatment decisions.

Continuous Glucose Monitoring System PDA: When the Continuous Glucose Monitoring System PDA is used together with compatible devices, it is intended to continuously measure the glucose in the interstitial fluid and is designed to replace fingerstick blood glucose (BG) testing for treatment decisions.

Indications

1. Type 1 and Type 2 Diabetes Mellitus
2. Special types of diabetes (excluding monogenic diabetes syndromes, diseases of the exocrine pancreas, and drug or chemical induced diabetes)
3. Abnormal blood glucose levels
4. Patients requiring improved glycemic control
5. People requiring frequent or continuous monitoring of blood glucose

Patients

Patients with diabetes (≥2 years old).

Intended User

The target users of this medical device are individuals aged 18 and above, who possess basic cognitive, literacy, and independent mobility

skills. It is intended for both medical professionals and non-professional adults who need to continuously or periodically monitor their own or others' glucose levels.

EN

Contraindications

The Continuous Glucose Monitoring System must be removed prior to Magnetic Resonance Imaging (MRI).
Don't wear your CGM sensor for computed tomography (CT) scan, or high-frequency electrical heat (diathermy) treatment.
Taking higher than the maximum dose of acetaminophen (e.g. > 1 gram every 6 hours in adults) may affect the rtCGMS readings and make them look higher than they really are.

The CGM System was not evaluated for the following persons:

- Pregnant women
- Peritoneal dialysis patients
- Patients with implanted pacemakers
- Patients with coagulation disorders or those taking anticoagulant drugs.

Warning

- Don't wear your CGM sensor for computed tomography (CT) scan, or high-frequency electrical heat (diathermy) treatment.

- Don’t wear your CGM while using electrocautery, electrosurgical units and diathermy equipment.
- The CGM system was not evaluated for the peritoneal dialysis patients, patients with implanted pacemakers and patients with coagulation disorders or those taking anticoagulant drugs. Before you use the Wellion SENSOR CGM system, review all the product instructions.
- The rtCGMS should not be used by Patients who have diffuse subcutaneous nodules.
- Before you use the Wellion SENSOR CGM system, review all the product instructions.
- The user’s manual includes all safety information and instructions for use.
- Talk to your healthcare professional about how you should use your sensor glucose information to help manage your diabetes.
- Failure to use the system according to the instructions for use may result in you missing a severe low blood glucose or high blood glucose event and/or making a treatment decision that may result in injury. If your glucose cautions and readings from the system do not match symptoms or expectations, use a fingerstick blood glucose value from a blood glucose meter to make diabetes treatment decisions. Seek medical attention when Appropriate.
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment

should be observed to verify that they are operating normally.

- Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Transmitter and PDA, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- After restarting your phone, please check again if Bluetooth is turned on. If it's turned off, please enable Bluetooth again to ensure real-time data transmission and notifications.
- Avoid these areas for sensor insertion:
 1. With loose skin or without enough fat to avoid muscles and bones.
 2. That get bumped or pushed, or that you lie on while sleeping.
 3. Within 3 inches of infusion or injection site.
 4. Near waistband or with irritations, scarring, tattoos, or lots of hair.
 5. With moles or scars.
- Android users, after enabling airplane mode, please double-check if Bluetooth is turned on. If it's turned off, please enable Bluetooth again to ensure real-time data transmission and notifications. iOS users don't need to consider this for the time being.

Precautions

- rtCGMS readings should only be used as a reference for the supplemental monitoring of diabetes mellitus and should not be used as a basis for clinical diagnosis.
- The rtCGMS should be completely removed before Magnetic Resonance Imaging (MRI).
- The rtCGMS contains many small parts that can be dangerous if swallowed.
- During rapid changes in blood glucose (more than 1.8 mg/dL or 0.1mmol/L per minute), glucose levels measured in interstitial fluid by the rtCGMS may not be the same as blood glucose levels. When blood glucose levels drop rapidly, the sensor may produce a higher reading than the blood glucose level; Conversely, when blood glucose levels rise rapidly, the sensor may produce a lower reading than the blood glucose level. When experiencing symptoms of hypoglycemia or hyperglycemia, the sensor readings are checked by using a blood glucose meter for a fingertip blood glucose test.
- Severe dehydration or excessive loss of water may result in inaccurate results. When you suspect you are dehydrated, consult a healthcare professional immediately.
- If you think the rtCGMS sensor reading is inaccurate or inconsistent, use a blood glucose meter to test your blood glucose level or calibrate the glucose sensor. If the problem persists, remove and replace the sensor.
- While extensive user testing was done on Wellion SENSOR CGMS in

Type I and Type II diabetics, the study groups did not include women with gestational diabetes.The performance of the rtCGMS has not been evaluated when used with another implantable medical device, such as a pacemaker.

EN

- Only Wellion SENSOR consumables should be used with the rtCGMS.
- If the product is not working properly or has been damaged, stop using the product.
- When used in a medical facility, the operator should wear gloves while using.

Potential clinical side-effects

Like any medical device, the Wellion SENSOR CGMS has potential side effects. The most common side effects include skin redness and skin ulceration at the sensor insertion site.

Additional security information



- Physiological differences between interstitial fluid and capillary whole blood may cause different glucose readings. Differences between sensor glucose readings from interstitial fluid and capillary blood can be observed during periods of rapid changes in blood glucose levels, such as after eating, insulin doses, or exercise.
- If you have to undergo a physical examination involving strong magnetic or electromagnetic radiation (e.g. MRI or CT), remove your sensor and insert a new sensor after the examination.
 The impact of these procedures on sensor performance has not been evaluated.



- The sensor Applicator is sterile in unopened and undamaged packages.
- Don't freeze the sensor. Do not use it after it expires.
- You are responsible for properly securing and managing your phone. If you suspect an adverse cyber security event related to the Wellion AiDEX SENSOR App, contact the customer service.
- Make sure that your phone and sensor kit are kept in a safe place, under your control. This is important to prevent anyone from accessing or tampering the System.
- The Wellion AiDEX SENSOR App is not intended for use on a phone that has been altered or customized to remove, replace or circumvent the manufacturer's Approved configuration or use restriction, or that otherwise violates the manufacturer's warranty.

Product components

Name: Wellion SENSOR continuous glucose monitoring system

System configuration: This product includes a Personal Diabetes Assistant, a Transmitter and a Glucose Sensor.
The Personal Diabetes Assistant will be referred to as PDA in the text below.

What you see	What it´s called	Model Number	What it does
 Glucose Sensor before insertion (Sensor Applicator)	Continuous glucose monitoring system sensor	G7-S01A: 10 days	• The Sensor-Applicator helps you to insert the Sensor under your skin. It contains a needle with which the skin is pierced and the flexible sensor is inserted into the skin. Once the sensor is in place, the needle is retracted into the container. • The Sensor is an Applied part which is only visible after Application, the sensor measures and stores glucose readings when worn on your body.
 Glucose Sensor after insertion			

What you see	What it´s called	Model Number	What it does
	Continuous glucose monitoring system transmitter	G7-T01A: 10 days	When the transmitter is used together with the glucose sensor base, when the sensor and transmitter are activated on one's body, it measures and stores glucose readings as well as transmits the readings to the PDA /the Wellion AiDEX Sensor App.
	Continuous glucose monitoring system PDA	G7-P02	The PDA is used to display glucose readings and historical data, when paired with an operational sensor-transmitter set.

Component description


Portable Diabetes Assistant

EN

The PDA may be purchased separately – information about its usage, refer to the separate PDA manual.

Wellion SENSOR Transmitter


The transmitter is used together with the glucose sensor base. When the sensor and transmitter are activated on one's body, it measures and stores glucose readings as well as transmits the readings to the PDA/the Wellion AiDEX SENSOR App.

Wellion SENSOR CGM package


The glucose sensor package contains a sterilized glucose sensor attached to a sensor base, and a sensor Applicator.

The package is sterilized by irradiation.

Using a new glucose sensor


Applying a sensor


- Recommended areas for sensor Application includes the abdomen, the outside and the back of the upper arm. For best performance, avoid excessive motion which may weaken the sensor and its adhesive tape. Avoid accidental knocking off the sensor. It is strongly recommended to put the sensor on your abdomen rather than arms or other areas.


EN

- Before Application, use alcohol pads to disinfect your selected skin area and wait for a minute to let the skin dry. Change your sensor sites regularly to avoid discomfort or skin irritation.
- Open the sensor package.



Note: Ensure the sensor is not expired by checking its expiration date. If the sensor is expired or the sterilized package is damaged, please do not use it

- Take out the sensor Applicator by turning it anticlock wise.
- Place the Applicator on the top of the desired sensor site. Press it firmly against skin and press the button to launch the sensor. After inserting the sensor, wait a few seconds to allow the sensor base patch to adhere to the skin.


- Remove the Applicator. Now, the sensor should be successfully attached.



EN


Note: Applying a sensor may cause bruising or bleeding. If bleeding does not stop, remove the sensor and Apply a new sensor to a different location.


- Align the transmitter to the sensor base and press the transmitter firmly onto the sensor base until it snaps onto the sensor base and attaches to it securely.



Note: If the orientation of the transmitter and sensor base mismatch, the transmitter cannot be installed properly.

Note: When the Sensor and Transmitter are combined, the protection level can reach to IP48, which means the level of protection against ingress of solid foreign objects is 4, and the level of protection against ingress of water with harmful effects is 8.

EN

- The non-recyclable Applicator is meant for single use only. Dispose the used sensor package and sensor Applicator according to your local regulations.
- Follow the steps described above carefully and be sure to only use components that are made by MicroTech Medical. Use of unapproved components can result in injury.

Pairing with a new transmitter

For pairing (and unpairing) your transmitter with the Wellion AiDEX Sensor App on your mobile phone, please refer to the App manual in the Help Center of the Wellion AiDEX Sensor App).

Sensor start-up

After the sensor is successfully Applied to your body and connected to the transmitter, the paired PDA/App will ask for a confirmation, if it is a new or a used sensor. If not automatically detected by the App, confirm a “new sensor”. Now the PDA/App will display “Sensor is warming up” for around 1 hour. After 1 hour you receive glucose values and alarms.

When the transmitter loses the connection with the sensor and reconnects, the App will ask you to confirm, if it is a new or an used sensor . If not automatically detected by the App, confirm a “used sensor”. The PDA/App will skip the warm-up process and starts displaying your current glucose value immediately.

Note: Do not cancel the warm-up process when using a new sensor. Skipping the new sensor initialization may make the system display incorrect glucose values.

Sensor calibration

Calibrating your rtCGMS using your BG meter is optional. You can manually input the BG meter value for calibration when the percentage difference is greater than 20%. It’s recommended to use BG meter readings of more than 4 mmol/L (72 mg/dL) to calibrate your rtCGMS .

Calibration with your BG meter may align your rtCGMS readings with your BG meter values.

For details about how to calibrate your rtCGMS, refer to the App manual in the Help Center of the Wellion AiDEX Sensor App.
If you want to calibrate your CGMS using your BG meter, enter the BG meter value within 5 minutes of taking a fingerstick.

Consistent BG meter: If you calibrate rtCGMS using your BG meter, use the same meter throughout the sensor session.


Only use fingertips test to calibrate rtCGMS. DO NOT use testing values from the alternative site (non-fingertip).

Calibration is not recommended when glucose values are rising or falling fast.

Important: Do not calibrate more than once every 6 hours.

Removing a sensor

- Remove the old sensor when the PDA/Wellion AiDEX SENSOR App indicates that the sensor has expired or if you feel any site irritation or discomfort.
- Carefully pull the tape that holds the sensor to the skin from the edge and slowly uncover it until the entire sensor is removed. You can use warm soapy water to remove any remaining sticky residue.
- Press the two locking arms on one side of the sensor base, pull the transmitter away from the sensor base, and then remove the transmitter so that it can be used again.
- The sensor is meant for one-time use only. Dispose the used sensor according to your local regulations.

Note: Please make sure the Wellion AiDEX SENSOR App is paired with the sensor before unpairing. If the sensor is not connected to the App, you may delete the sensor record permanently by clicking "Delete".

Replacing the sensor

After sensor expiration, your sensor will automatically stop working and needs to be replaced. In addition, if you notice any irritation or discomfort at the Application site, or if the Application fails, you should replace your sensor.

Wellion AiDEX SENSOR App – Diabetes Management System

The Wellion AiDEX SENSOR App is used in combination with the Wellion SENSOR CGMS for continuous glucose monitoring.
For more information visit www.wellion.at

Caring for your Continuous Glucose Monitoring System (rtCGMS)

Cleaning the transmitter

- Disinfect the transmitter with an alcohol wipe.
- Do not use solvents, nail polish remover, or paint thinner to scrub the outer surface.
- Keep the transmitter dry.

- Do not use any lubricant.

Disposal

EN

Dispose transmitters and sensors that you do not use any more in compliance with local regulations for electronic devices, batteries, sharps, and biohazard materials. Please do not discard old products or accessories directly. For further information on how to dispose of system components, please contact customer service.

Transportation

Note: Do not immerse the transmitter in water or other liquids.

Avoid placing heavy weight on top of the transmitter. Avoid direct sunlight and rain.

Storage

If you are temporarily not using the transmitter, or sensor system, store the components in a cool, dry, clean and well-ventilated area.
The transmitter and sensor are precision instruments. If they fail, they can only be returned to the manufacturer to repair. No third-party individuals or organizations are allowed to perform repairs. Circuit diagrams and component lists are not provided in the manual.

Specifications

<table>
<tr><th colspan="2" rowspan="2">Item</th><th colspan="3">Subcomponent</th></tr>
<tr><th>Transmitter</th><th>Sensor</th><th>PDA</th></tr>
<tr><td rowspan="3">10 days</td><td rowspan="3">Model number
&
Use Life</td><td>WELL19-S10D00T
(G7-T01A)</td><td>WELL19-S10D
(G7-S01A)</td><td>WELL19-S00PDA
WELL19-S00PDAMM
(G7-P02)</td></tr>
<tr><td rowspan="2">4 years</td><td>10 days</td><td rowspan="2">4 years</td></tr>
<tr><td>Shelf life: 1 year</td></tr>
<tr><td colspan="2">Operating temperature</td><td colspan="3">5-40°C (41-104°F)</td></tr>
<tr><td colspan="2">Operating humidity</td><td colspan="3">10-93% (non-condensing)</td></tr>
<tr><td colspan="2">Storage and transportation temperature</td><td>-20 - 60°C</td><td>4 - 30°C</td><td>-20 - 60°C</td></tr>
<tr><td colspan="2">Storage and transportation humidity</td><td colspan="3">5-95% (non-condensing)</td></tr>
<tr><td colspan="2">Storage and transportation pressure</td><td colspan="3">700hpa-1060hpa</td></tr>
<tr><td colspan="2">Ingress protection level</td><td colspan="2">IP48</td><td>IP20 (not waterproof)</td></tr>
<tr><td colspan="2">Detection range</td><td colspan="3">2,0 mmol/L-25,0 mmol/L (36 mg/dL - 450 mg/dL)</td></tr>
<tr><td colspan="2">Wireless frequency and bandwidth</td><td colspan="3">Frequency: 2,402GHz-2,48 GHz
Bandwidth: 1Mbps</td></tr>
<tr><td colspan="2">Wireless modulation</td><td colspan="3">GFSK</td></tr>
<tr><td colspan="2">Radiated power</td><td colspan="3">-2dBm</td></tr>
</table>

EN

Electromagnetic compatibility

EN

These devices are intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that the device is used in such an environment.

Portable and mobile RF communication interference may have an impact on the device.

The use of accessories other than those specified for the device is not recommended. They may result in increased emissions or decreased immunity of the device.

The device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device should be observed to verify normal operation in the configuration in which it will be used.

The basic performance is described in the following table:

Performance	Description
Measurement accuracy	When the glucose concentration is >4.2 mmol/L (75mg/dL), the sensor accuracy deviation does not exceed ±20%; When the glucose concentration is ≤4.2 mmol/L (75mg/dL), the accuracy deviation does not exceed ± 1,1 mmol/l (20mg/dL)

IEC 60601-1-2:

Guidance and manufacturer's declaration – electromagnetic immunity		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should ensure that it is used in such an environment.		
Emissions test	**Compliance**	**Electromagnetic environment - guidance**
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply.
Harmonic emissions IEC 61000-3-2	Class A	Move to a place within the normal operating temperature range and repeat the test.
Voltage fluctuations/ Flicker emissions IEC 61000-3-3	Complies	Repeat test. If you see the same result, contact your healthcare professional immediately.

IEC 60601-1-2:

Guidance and manufacturer's declaration – electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should ensure that it is used in such an environment.			
Immunity test	**IEC 60601 test level**	**Compliance level**	**Electromagnetic environment - guidance**
Electrostatic discharge(ESD) IEC 60601-4-2	± 6 kV Contact ± 8 kV Air	± 6 kV Contact ± 8 kV Air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient burst IEC 61000-4-4	± 2 kV power cord ± 1 kV input/output	± 2 kV power cord ± 1 kV input/output	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1KV line to GND ±2KV line to GND	±1KV line to GND ±2KV line to GND	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4- 11	<5% UT for 0.5 weeks (>95% dip in UT) 40% UT for 5 weeks (60% dip in UT) 70% UT for 25 weeks (30% dip in UT) <5% UT for 5s (>95% dip in UT)	<5% UT for 0.5 weeks (>95% dip in UT) 40% UT for 5 weeks (60% dip in UT) 70% UT for 25 weeks (30% dip in UT) <5% UT for 5s (>95% dip in UT)	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	The power frequency magnetic field should have the characteristics of power frequency magnetic field level at a atypical place in a typical commercial and hospital environment.

IEC 60601-1-2:

Guidance and manufacturer's declaration – electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should ensure that it is used in such an environment.			
Immunity test	**IEC 60601 test level**	**Compliance level**	**Electromagnetic environment - guidance**
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3V(Vrms) 150 kHz ~ 80 MHz 10V(engi-neering medical frequency band) 150kHz~80MHz 3 V/m 80MHz-2,7GHz	3V(Vrms) 10V (engineering medical frequency band) 3V/m	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation Applicable to the frequency of the transmitter. Recommended separation distance $d=1.2\sqrt{P}$ $d=1.2\sqrt{P}$ 80MHz~800MHz $d=2.3\sqrt{P}$ 800MHz~2.5GHz where P is the maximum output power rating of the transmitter in watts (W) according to the trans-mitter manu- facturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 

EN

IEC60601-1-2:

EN

Recommended separation distances between portable and mobile RF communications equipment and the device			
These devices are intended for use in an environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment(transmitters) and the device as recommended below, according to the maximum output power of the communications equipment			
Maximum rated output power of transmitter (W)	**Separation distance according to frequency of transmitter (m)**		
	150kHz-80MHz $d=1.2\sqrt{P}$	80MHz-800MHz $d=1.2\sqrt{P}$	800MHz-2.5GHz $d=2.3\sqrt{P}$
0,01 0,1 1 10 100	0,12 0,38 1,2 3,8 12	0,12 0,38 1,2 3,8 12	0,23 0,73 2,3 7,3 23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation Applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range Applies. NOTE 2 These guidelines may not Apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Clinical performance

To demonstrate the accuracy performance of the CGM System, two prospective clinical studies were conducted at 9 centers in China. The studies included both adult (18 years and older) and pediatric (2 to 18 years) participants. The studies enrolled a total of 121 adult (18 years and older) and 90 pediatric (2 to 18 years) participants with Type 1 diabetes mellitus or Type 2 diabetes mellitus. A total of 209 subjects completed the studies.

Participants wore either two or four sensors for up to 14 days. At least 2 sensors were Applied to a subset of participants for the precision study to compare the variability of readings between sensors. Adult participants wore their rtCGMS in the abdomen and upper arm; pediatric subjects wore their rtCGMS in the abdomen. Clinic session(s) took place on Day 1-2, Day 6-8 and Day 14-15. Depending on the participant's age, they participated in one clinic session of varying duration.

Accuracy of the rtCGMS was evaluated by comparing:

rtCGMS values obtained at the same or similar time point to the comparator device YS1 in the adult study.

rtCGMS values obtained at the same or similar time point to the comparator device EKF in the pediatric study.

Accuracy results obtained from Adult Study (Home user)

Indicator	Result
Mean Absolute Relative Difference (MARD%)	8.66%
When blood glucose level > 4.4 mmol/L (79 mg/dL)	
Results within a deviation range of ±20% from the reference value.	94.56%
When blood glucose level ≤ 4.4 mmol/L (79 mg/dL)	
Results within a deviation range of ±1.1mmol/L (19mg/dL) from the reference value.	1.72%
The percentage of data points that fall within Clarke error grid zones A+B	100%
The percentage of data points that fall within Consensus error grid zones A+B	100%

Caution rate

Success rate of hyperglycemic caution: 89.23% (hyperglycemic caution threshold set at 11.1mmol/L (200mg/dL));
Success rate of hypoglycemic caution: 100.0% (hypoglycemic caution threshold set at 4.4mmol/L (79mg/dL)).

Accuracy results obtained from Adult Study (Hospital user)

Indicator	Result
Mean Absolute Relative Difference (MARD%)	9.55%
When blood glucose level > 4.4 mmol/L (79 mg/dL)	
Results within a deviation range of ±20% from the reference value.	93.02%
When blood glucose level ≤ 4.4 mmol/L (79 mg/dL)	
Results within a deviation range of ±1.1mmol/L (19mg/dL) from the reference value.	2.44%
The percentage of data points that fall within Clarke error grid zones A+B	100%
The percentage of data points that fall within Consensus error grid zones A+B	100%

Caution rate

Success rate of hyperglycemic caution: 90.57% (hyperglycemic caution threshold set at 11.1mmol/L (200mg/dL));
Success rate of hypoglycemic caution: 97.37% (hypoglycemic caution threshold set at 4.4mmol/L (79mg/dL)).

Safety

3 device-related adverse events (bleeding, pain and skin redness at the insertion site) were recorded in both user groups.

EN

Accuracy results obtained from pediatric Study

Indicator	Result
Mean Absolute Relative Difference (MARD%)	9.51%
When glucose concentration ≥ 3.90mmol/L (70 mg/dL) and< 10.00mmol/L (180 mg/dL)	
Results within a deviation range of ± 15% from the reference value.	82.10%
Results within a deviation range of ± 40% from the reference value.	99.93%
When the glucose concentration ≥ 10.00mmol/L (180 mg/dL)	
Results within a deviation range of ± 15% from the reference value.	89.48%
Results within a deviation range of ± 40% from the reference value.	100.00%
When the glucose concentration < 3.90mmol/L (70 mg/dL)	
Results within a deviation range of ± 0.83mmol/L (14 mg/dL) from the reference value.	91.28%
Results within a deviation range of ± 2.22mmol/L (40 mg/dL) from the reference value.	91.28%
The percentage of data points that fall within Clarke error grid zones A+B	98.85%
The percentage of data points that fall within Consensus error grid zones A+B	99.95%

Caution rate

Success rate of hyperglycemic caution: 90.57% (hyperglycemic caution threshold set at 11.1mmol/L (200mg/dL));
Success rate of hypoglycemic caution: 97.37% (hypoglycemic caution threshold set at 4.4mmol/L (79mg/dL)).




Clinical side effect

6 device-related adverse events (pain, itchy skin, swelling, pus thrombus, subcutaneous mass, and dermatitis at the sensor insertion site) were recorded in the pediatric trial.

Appendix

Symbols

EN

	Do not re-use
	Type BF Applied part
	Temperature limit
	Atmospheric pressure limitation
	Humidity limitation
	Single sterile barrier system with protective packaging outside using irradiation
	Refer to instruction manual
	Medical device
	Non-ionizing radiation
	Caution
	Consult the Instructions for Use
	Manufacturer

Importer

EU REP

Authorised Representative in the European Community

MR unsafe

Do not use if package is broken

Date of manufacture

Use-by date

LOT

Batch Number

SN

Serial Number

Waste Electrical and Electronic Equipment (WEEE)

UDI

Unique device identifier

CE 0197

CE Mark

IP48

The level of protection against ingress of solid foreign objects is 4 (Protected against solid foreign objects of 1,0 mm Ø and greater).
The level of protection against ingress of water with harmful effects is 8 (Protected against the effects of continuous immersion in water).

Potential interference information

EN

It has been studied that when users take normal doses of ascorbic acid or acetaminophen (ascorbic acid blood concentration <0.055mmol/L (1mg/dL), acetaminophen blood concentration <0.055mmol/L (1mg/dL)); the drug will not interfere with the sensor glucose measurement. The hydroxyurea has a significant impact on CGM measurement values. The error size depends on the actual concentration of the blood uric acid value. If the user feels that the current physical condition does not match the glucose readings obtained by the Continuous Glucose Monitoring System or suspects that the measurements may be inaccurate, the blood glucose test can be performed using a finger blood glucose meter and corresponding management actions can be taken based on the test values. When using a finger blood glucose meter, record your blood glucose values promptly after measurement to avoid forgetting or inaccuracies in the readings.

Potential Risks

Inaccurate glucose values

Exposure to heat for longtime may cause in accurate results.

Mild to severe to sensor related-wear reactions

E.g. allergic reaction, moderate to severe itching, rash, erythema,

bleeding, minor infection at the insertion site, discomfort during insertion.

Hyperglycemia or hypoglycemia

Hypo- and hyperglycemia events stemming from missed cautions or sensor inaccuracies.

Potential clinical benefits

The potential clinical benefits of the Wellion SENSOR CGM system are improved management of A1C and TIR for tighter glycemic control. Shortened time spent in hypoglycemia and hyperglycemia reduction in hypo- and hyperglycemia events in diabetes patients.

Glossary

EN

Blood glucose meter

A device used to measure the levels of glucose in the blood.

Blood glucose results

The concentration of glucose in the blood is measured as either milligrams of glucose per deciliter of blood (mg/dL) or millimol of glucose per liter of blood (mmol/L).

Continuous glucose monitor (CGM)

A CGM uses a small sensor inserted below your skin to measure the amount of glucose in the fluid in your skin, called interstitial fluid. Those glucose results are then sent to an PDA/App, where they are displayed as glucose levels and long-term glucose trends.

Hyperglycemia (high blood glucose)

High levels of glucose in the blood, also known as high blood glucose. When left untreated, hyperglycemia can lead to serious complications. Talk to your healthcare professional to determine your high glucose level.

Hypoglycemia (low blood glucose)

Low levels of glucose in the blood, also known as low blood glucose. When left untreated, hypoglycemia can lead to serious complications.

Talk to your healthcare professional to determine your low glucose level.

Interstitial fluid
The fluid that surrounds all the cells of the body.

Insulin
A hormone produced by the pancreas that regulates the metabolism of glucose and other nutrients. If the pancreas is damaged and does not produce insulin, a doctor may prescribe insulin injections to help people with diabetes process glucose (sugar).

Limitations
A safety statement outlining specific situations in which the Wellion SENSOR CGM should not be used because it may be harmful to you or damage the system.

mg/dL
Milligrams per deciliter; one of two standard units of measurement for the concentration of blood glucose (sugar).

mmol/L
Millimol per liter; one of two standard units of measurement for the concentration of blood glucose (sugar).

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MicroTech Medical (Hangzhou) Co., Ltd.