wellion

Comtest 4in1 REF WELLCOMBO

(DE) Wellion Combotest 4in1 (SARS-CoV-2/Influenza A+B/RSV Ag Combo Schnelltest) Packungsbeilage

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INTENDED USE
The Wellion Combotest 4in1 (SARS-CoV-2/Influenza A+B/RSV Ag Combo Rapid Test) is a rapid chromatographic immunoassay for the qualitative detection of SARS-CoV-2, Influenza A, Influenza B and Respiratory Syncytial Virus (RSV) antigens in nasopharyngeal swab specimens from persons with suspected SARS-CoV-2/Influenza/RSV infection in conjunction with clinical presentation and other laboratory test results.
The test is for in vitro diagnostic use only. For professional use only. It is intended for clinical laboratories and healthcare professional use only for point-of-care testing. Not for at-home testing.

PRINCIPLE
The SARS-CoV-2 Rapid Test has one line of anti-SARS-CoV-2 antibody on the detection line (T line) and one line of anti-mouse IgG antibody on the quality control line (C line). When extracted specimen is added to the specimen well it will react with the labeled antibody to form a complex; the mixture then migrates through the membrane by capillary action and interacts with the anti-SARS-CoV-2 antibody on the detection line. If the specimen contains SARS-CoV-2 antigen, the detection line will appear, indicating the presence of SARS-CoV-2 antigen. Otherwise, the test result will be negative. The test device also contains a quality control line C, which should appear for all valid tests. If the quality control line C does not appear, the test result will be invalid even if the detection line appears.
The Influenza A+B Rapid Test has one line of anti-influenza A antibody on the detection line (A line), one line of anti influenza B antibody on the detection line (B line) and one line of anti-mouse IgG antibody on the quality control line (C line). When the specimen is added to the specimen well, it will react with the labeled antibody to form a complex; the mixture then migrates through the membrane by capillary action and interacts with the anti-influenza A antibody and anti-influenza B antibodies on the detection line. If the specimen contains influenza A or influenza B antigens, the detection line will appear and indicating the presence of influenza A or influenza B antigens. Otherwise, the test result will be negative. The test device also contains a quality control line C, which should appear for all valid tests. If the quality control line C does not appear, the test result will be invalid even if the detection line appears.
The RSV Rapid Test has one line of anti-RSV antibody on the detection line (T line) and one line of anti-mouse IgG antibody on the quality control line (C line). When extracted specimen is added to the specimen well, it will react with the labeled antibody to form a complex; the mixture then migrates through the membrane by capillary action and interacts with the anti-RSV antibody on the detection line. If the specimen contains RSV antigen, the detection line will appear red indicating the RSV antigen is positive. Otherwise, the test result will be negative. The test device also contains a quality control line C, which should appear for all valid tests. If the quality control line C does not appear, the test result will be invalid even if the detection line appears.

COMPOSITION
Each test kit contains test, sterile swab, extraction solution (in the sealed tube), tube holder, sterile swab and package insert.
Materials required but not provided: timer.

STORAGE AND STABILITY
Store the test kit at a cool, dry place between 36-86°F (2-30°C). Keep away from light. Exposure to temperature and / or humidity outside the specified conditions may cause inaccurate results.
Do not freeze. Use the test kit at temperatures between 59-86°F (15-30°C).
Use the test kit between 10-90% humidity.
Do not use the test kit beyond the expiration date (printed on the foil pouch and box).

Note: All expiration dates are printed in format YYYY-MM-DD. 2022-06-18 indicates June 18, 2022.2022-06-18.

WARNINGS AND PRECAUTIONS
- For in vitro diagnostic use only.
- This package insert must be read completely before performing the test. Failure to follow the instructions may cause inaccurate results.
- Do not open the sealed pouch unless ready to conduct the assay.
- Do not use expired devices.
- Bring all reagents to room temperature (15°C-30°C) before use.
- Do not use the components in any other type of test kit as a substitute for the components in this kit.
- Wear protective clothing and disposable gloves while handling the kit reagents and clinical specimens. Wash hands thoroughly after performing the test.
- Do not smoke, drink, or eat in areas where specimens or kit reagents are being handled.
- Dispose of all specimens and materials used to perform the test as biohazardous waste.
- Handle the Negative and Positive Control in the same manner as patient specimens.

LIMITATIONS
- This test is for in -vitro diagnostic use only and cannot be re-used.
- A positive test should be treated as potentially infectious materials and should be disposed of properly.
- The test kit should be protected away from direct sunlight, moisture and heat.
- Please check if the test kit has any damage and check the expiry date before use.
- Adequate volume will affect the accuracy of the test result. Inaccurate sample volume may cause a false positive or negative result.
- The test result must be evaluated in conjunction with other clinical data available to the physician.
- Positive test results do not rule out co-infections with other pathogens.
- Positive and negative predictive values are highly dependent on prevalence. The local prevalence should be taken into consideration when interpreting diagnostic test results.
- Please be very careful when collecting nasopharyngeal swab specimen from children.
- Components from different batches are not allowed to be used in combination.

SPECIMEN COLLECTION
Consider any materials of human origin as infectious and handle them according established procedures.

Specimen Collection
Use the nasopharyngeal swab supplied in the kit. It is important to obtain as much secretion as possible. Therefore, to collect a nasopharyngeal swab sample, carefully insert the sterile swab into the nostril that presents the most secretion under visual inspection. Keep the swab near the nasal septum floor of the nose while gently pushing the swab into the posterior nasopharynx. Rotate the swab several times, and then remove it from the nasopharynx.

TEST PROCEDURE
Allow the Test Devices and Extraction Solution to equilibrate to 15-30°C. Bring to room prior to use.

1. Open the extraction solution (in the sealed tube).

2. Collect specimen refer to Specimen Collection.
3. Insert the swab with collected specimen into the extraction tube filled with the extraction solution. Roll the swab 5 times while pressing the head against the bottom and side of the extraction tube.
4. Remove the swab while squeezing the sides of the tube to extract the liquid from the swab. Try to release as much liquid as possible. Dispose the used swab as a biohazard waste.

5. Put on the tube tip.

6. Take out a test device from sealed foil pouch and put it on a clean and level surface.
7. Apply 3 drops of the extracted specimen into the each specimen well. Please avoid bubbles during applying.

8. Read the test result at 15 minutes. Don't read the result after 20 minutes.

INTERPRETATION OF TEST RESULTS
For SARS-CoV-2:

Positive Results:
Both the quality control line C and the detection line T appear.
Negative Results:
Only the quality control line C appears, with no other line appearing on the detection area.
Invalid Result:
Quality control line C fails to appear indicating the test is invalid, no matter if the detection line appears or not. Collect a new specimen and perform another test with a new test device.

For Influenza A+B:

Positive Results:
Positive influenza A antigen:
Both the quality control line C and the influenza A detection line appear, while the influenza B detection line does not appear.
Positive influenza B antigen:
Both the quality control line C and the influenza B detection line appear, while the influenza A detection line does not appear.
Positive influenza A and B antigen:
All 3 lines appear, including the quality control line C and the influenza A and influenza B detection lines.
Negative Results:
Only the quality control line C appears, with no other line appearing on the influenza A and influenza B detection lines. It indicates the test result is negative for both influenza A and influenza B antigens.
Invalid Results:
Quality control line C fails to appear indicating the test is invalid, no matter if the influenza A or influenza B detection lines appear or not. Collect a new specimen and perform another test with a new test device.

For RSV:

Positive Results:
Both the quality control line C and the detection line T appear.
Negative Results:
Only the quality control line C appears, with no other line appearing on the detection area.
Invalid Results:
Quality control line C fails to appear indicating the test is invalid, no matter if the detection line appears or not. Collect a new specimen and perform another test with a new test device.

QUALITY CONTROL
Internal procedural controls are included in the test. A colored line appearing in the control region (C) is the internal procedural control line. The procedural control line indicates that sufficient flow has occurred, and the functional integrity of the test device has been maintained. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested

PERFORMANCE
Clinical Sensitivity/Clinical Specificity
A total of 562 specimens were tested using the Wellion Combotest 4in1 (SARS-CoV-2/Influenza A+B/RSV Ag Combo Rapid Test). These specimens were obtained by nasopharyngeal swab from symptomatic patients. The performance of the Wellion Combotest 4in1 (SARS-CoV-2/Influenza A+B/RSV Ag Combo Rapid Testwas compared to a commercialized molecular assay.

Table 1: SARS-CoV-2 Performance of the Wellion Combotest 4in1 (SARS-CoV-2/Influenza A+B/RSV Ag Combo Rapid Test) as Compared to PCR Test

Wellion Combotest 4in1 (SARS-CoV-2/ Influenza A+B/RSV Ag Combo Rapid Test) - for SARS-CoV-2	PCR		
	Positive	Negative	Total
Positive	162	1	163
Negative	4	255	259
Total	166	256	422
Sensitivity	96,23% (158/166, 95%CI, 90,70%–98,52%)		
Specificity	99,61% (255/256, 95%CI, 97,82%–99,93%)		
Accuracy	98,67% (157/162, 95%CI, 96,81%–99,47%)		

Table 2: Influenza A Performance of the Wellion Combotest 4in1 (SARS-CoV-2/Influenza A+B/RSV Ag Combo Rapid Test) as Compared to PCR Test

Wellion Combotest 4in1 (SARS-CoV-2/ Influenza A+B/RSV Ag Combo Rapid Test) - for Influenza A	PCR		
	Positive	Negative	Total
Positive	55	4	59
Negative	3	500	503
Total	58	504	562
Sensitivity	94,83% (55/58, 95%CI, 85,86%–98,25%)		
Specificity	99,21% (500/504, 95%CI, 98,04%–99,66%)		
Accuracy	98,58% (556/562, 95%CI, 96,93%–99,29%)		

Table 3: Influenza B Performance of the Wellion Combotest 4in1 (SARS-CoV-2/Influenza A+B/RSV Ag Combo Rapid Test) as Compared to PCR Test

Wellion Combotest 4in1 (SARS-CoV-2/ Influenza A+B/RSV Ag Combo Rapid Test) - for Influenza B	PCR		
	Positive	Negative	Total
Positive	32	1	33
Negative	1	528	529
Total	33	529	562
Sensitivity	96,97% (32/33, 95%CI, 84,68%–99,46%)		
Specificity	99,81% (528/529, 95%CI, 98,96%–99,98%)		
Accuracy	99,64% (560/562, 95%CI, 98,72%–99,96%)		

Table 4: RSV Performance of the Wellion Combotest 4in1 (SARS-CoV-2/Influenza A+B/RSV Ag Combo Rapid Test) as Compared to PCR Test

Wellion Combotest 4in1 (SARS-CoV-2/ Influenza A+B/RSV Ag Combo Rapid Test) - for RSV	PCR		
	Positive	Negative	Total
Positive	111	2	113
Negative	2	447	449
Total	113	449	562
Sensitivity	98,23% (111/113, 95%CI, 93,28%–99,53%)		
Specificity	99,55% (447/449, 95%CI, 97,92%–99,93%)		
Accuracy	98,94% (556/562, 95%CI, 97,78%–99,57%)		

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