wellion®

STREP A Rapid Test

REF: VIST-702

EN **Wellion STREP A Rapid Test** / Instruction for Use

INTENDED PURPOSE
The Strep A Rapid Test is a lateral flow chromatographic immunoassay for the qualitative detection of group A Streptococcus antigens in human throat swab samples. It is intended to be for auxiliary diagnosis of group A hemolytic streptococcus infection. For *in vitro* diagnostic use only.

SUMMARY
Streptococcus pyogenes is a species of non - motile gram - positive cocci which are categorized as Lancefield group A, which may cause serious infections such as pharyngitis, respiratory infection, impetigo, endocarditis, meningitis, puerperal sepsis, and arthritis.[1] Left untreated, these infections can lead to serious complications, including rheumatic fever or a peritonsillar abscess.[2] Traditional identification procedures for group A streptococcal infection involve the isolation and identification of viable organisms using techniques that require 24 to 48 hours or longer.[3,4] The Strep A Rapid Test is an immunological test designed to specifically identify Streptococcus A bacteria in oropharyngeal infections, providing results within 5 minutes. The test utilizes antibodies specific for whole cell Lancefield Group A Streptococcus to selectively detect Strep A antigen in a throat swab specimen.

TEST PRINCIPLE
The Strep A Rapid Test is a lateral flow immunoassay based on the principle of the double-antibody sandwich technique. The test cassette consists of: 1) a burgundy conjugate pad containing rabbit anti-group A Streptococcus monoclonal antibodies conjugated with colloidal gold (Strep A Ab conjugates), 2) a nitrocellulose membrane strip containing a test line (T line) and a control line (C line). The T line is pre-coated with further rabbit anti-group A Streptococcus monoclonal antibodies, and the C line is pre-coated with goat anti-rabbit IgG.
When an adequate volume of a test specimen is dispensed into the sample well of the test cassette, the specimen migrates by capillary force along the membrane. The group A Streptococcus antigens, if present in the specimen, will bind to the Strep A Ab conjugates. The immunocomplex is then captured on the membrane by the pre-coated rabbit anti-group A Streptococcus monoclonal antibody reagent, forming a burgundy T line, indicating a group A Streptococcus positive test result.
The test contains an internal control (C line) which, regardless of colour development on the T line, should exhibit a burgundy coloured control line, caused by the immunocomplex of the control antibodies. If no C line appears, the test result is invalid and the specimen must be retested using another test cassette.

WARNINGS AND PRECAUTIONS
- For in-vitro diagnostic use only.
- All specimens should be considered capable of transmitting disease. Take appropriate precautions when collecting, handling, storing, and disposing of patient samples and used kit contents. And follow biosafety guidelines.
- Wear appropriate personal protective equipment (eg gown, gloves, goggles) when handling the contents of this kit.
- Proper sample collection, storage, and transportation are critical to the performance of this test.
- Dispose after use. The test cassette is not suitable for reuse.
- Reagent 2 contains an acidic solution. If the solution comes into contact with skin or eyes, flush with plenty of water.
- Do not touch the reaction area on the test cassette.
- Do not use beyond the expiry date.
- Do not use the test cassette if the pouch is punctured or poorly sealed.
- Do not interchange bottle caps between reagents.
- Do not interchange external control solution caps.
- Do not smoke, eat or drink in the testing area.
- Testing should be performed by professionally trained personnel working in an accredited laboratory or clinic and samples should be collected by qualified medical personnel.
- Test results should be interpreted by a physician along with clinical findings and other laboratory test results .
- Disposal of diagnostic reagents: All specimens and used kits carry a risk of infection. Disposal of diagnostic reagents must comply with local infectious disease disposal laws or laboratory regulations.


COMPONENTS
- The Strep A Rapid Test contains:
- 1x test cassette (packaged in pouch with desiccant)
- 1x tongue depressor
- 1x disposable swab
- 1x reagent 1 & 1x reagent 2
- 1x test tube & 1x tube tip
- 1x tube stand
- 1x Instructions for use

MATERIALS REQUIRED BUT NOT SUPPLIED
- Timer
- Personal protective equipment, such as protective gloves, medical masks, lab coats, etc.
- Appropriate biohazardous waste containers and disinfectants.

STORAGE AND STABILITY
- Store the test kit between 2-30°C. Keep away from light. Exposure to temperatures and / or humidity outside the specified conditions may cause inaccurate results.
- Do not freeze. Use the test at temperatures between 15-30°C.
- Do not use beyond the expiration date (printed on the foil pouch and box).


Note: All expiration dates are printed in Year-Month-Day format. 2022-06-18 indicates June 18, 2022.

SPECIMEN COLLECTION AND PREPARATION


1. **Specimen collection**
Collect a throat swab sample using the disposable swab provided with the test kit. The sampler should secure the tongue with a tongue depressor, then using the swab to pass beyond the root of the tongue to the posterior pharyngeal wall, tonsil crypts, and lateral walls. Gently wipe the throat area 3-5 times with the swab to collect mucosal cells, avoiding touching with the tongue, inner cheeks, and teeth. Inserting the swab back into the sampling device or an appropriate transport container. Alternatively, swabs containing modified Stuart or Amies medium can be used for sample collection.
2. **Specimen handling**
Testing is recommended immediately after collecting a throat swab sample. If testing cannot be performed immediately, throat swab samples should be stored in a clean, dry, airtight container at room temperature for up to 4 hours and at 2-8°C for up to 24 hours. DO NOT FREEZE.
If culture is required, throat swab samples should be transferred to Group A Selective (GAS) blood agar plates for culture.

TEST PROCEDURE
Allow the tests, samples and reagents to come to room temperature (15-30° C) before testing. Remove the test device from the sealed foil pouch and use it as soon as possible. Best results will be obtained if the test is performed immediately after opening the foil pouch.

1. Add 6 drops of Reagent 1 (white bottle cap) and 4 drops of Reagent 2 (red bottle cap) to the test tube and mix properly by gently swirling the test tube.
2. Collect specimen refer to Specimen Collection.
3. After collecting the sample, insert the swab's cotton tip into the plastic test tube placed in the tube stand, containing the previously added two reagents.
4. Hold the plastic handle of the swab firmly and rotate the swab about 10 times against the sides of the test tube to mix the solution thoroughly. Leave the swab to incubate for 1 minute.
5. Remove the test tube (with the swab and extraction fluid still in it) from the tube stand. Using your thumb and index finger, press the sides of the test tube to release as much fluid as possible from the swab's cotton tip and collect it in the test tube. Remove the swab.
6. Dispose of the swab according to local regulations and put the test tube back into one of the openings of the tube stand.
7. Add the tube tip to the plastic test tube.
8. Add 3 drops to the sample well of the test cassette.
Note: If the dispensed drop contains air bubbles, add another drop to the sample well.
9. Wait for the coloured line(s) to appear. Read the result after 5 minutes. Do not read the results after more than 10 minutes.

RESULT INTERPRETATION

Positive*: Two distinct coloured lines appear. One line should be in the control line region (C) and another line should be in the test line region (T).
*Note: The intensity of the colour in the test line region (T) may vary depending on the concentration of group A Streptococcus antigens present in the specimen. Therefore, any shade of colour in the test line region (T) should be considered positive.
Negative: A coloured line appears in the control line region (C), no line appears in the test line region (T), indicating that the result is negative.
Invalid: The control line fails to appear. Insufficient sample volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test cassette. If the problem persists, stop using the test cassette and contact your local distributor.

QUALITY CONTROL
Internal procedural controls are included in the test. A color line appearing in the control region (C) is an internal positive procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

LIMITATIONS
- The Strep A Rapid Test is intended for professional, in-vitro diagnostic use only and is intended for the qualitative detection of group A Streptococcus antigens in throat swab samples.
- A positive result of the Strep A Rapid Test only indicates the presence of group A Streptococcus antigens in the sample and should not be used as the sole criterion for the diagnosis of a group A hemolytic streptococcal infection.
- Negative results should be confirmed by culture. A negative result may be obtained if the concentration of group A Streptococcus antigens in the throat swab is insufficient or below the detection limit of the test.
- Excessive blood or mucus in the swab sample can impair the test performance and lead to false positive results.
- As with all diagnostic tests, a doctor can only make a definitive diagnosis after evaluating all clinical and laboratory findings.

PERFORMANCE CHARACTERISTICS
1. **Accuracy**
Based on 247 samples, the results of the Strep A Rapid Test compared to a culture method are shown as follows:

Strep A Rapid Test	Culture Method		
	Positive	Negative	Total
Positive	86	2	88
Negative	4	155	159
Total	90	157	247
Diagnostic sensitivity	95.56% (86/90, 95%CI, 89.12%~98.26%)		
Diagnostic specificity	98.73% (155/157, 95%CI, 95.47%~99.65%)		
Diagnostic accuracy	97.57% (241/247, 95%CI, 94.80%~98.88%)		

2. **Cross-reactivity**
There was no cross-reaction with potential cross-reactive substances. The following organisms were tested at 1 × 10⁷ organisms/swab and showed negative results.

Potential cross-reactive substances	Potential cross-reactive substances	Potential cross-reactive substances
Group F *Streptococcus*	Group B *Streptococcus*	*Streptococcus pneumoniae*
Streptococcus mutans	*Staphylococcus aureus*	*Corynebacterium diphtheriae*
Candida albicans	*Pseudomonas aeruginosa*	*Neisseria meningitidis*
Neisseria sicca	*Moraxella catarrhalis*	*Group C Streptococcus*
Group G Streptococcus	*Streptococcus sanguinis*	*Enterococcus faecalis*
Staphylococcus epidermidis	*Serratia marcescens*	*Klebsiella pneumoniae*
Bordetella pertussis	*Neisseria gonorrhoeae*	*Neisseria subflava*
Haemophilus influenzae	/	/

3. **Interfering Substances**
There was no interference with the potentially interfering substances listed below.

Interfering substance	Concentration in specimen	Interfering substance	Concentration in specimen
4-Acetamidophenol	10 mg/mL	Benzocaine throat spray (Cepacol)	5% by volume
Menthol Throat Lozenges	5% w/v	Blood, type A	2% (v/v)
Acetylsalicylic acid	20 mg/mL	Blood, type B	2% (v/v)
Albuterol	0.083 mg/mL	Blood, type AB	2% (v/v)
Amantadine	500 ng/mL	Blood, type O	2% (v/v)
Mometasone	500 ng/mL	Dexamethasone	10 mg/mL
Mouthwash Listerine	5% (v/v)	Dextromethorphan	10 mg/mL
Mouthwash Lion	5% (v/v)	Diphenhydramine HCl	5 mg/mL
Ascorbic acid chewable tablets	5% by weight	Ibuprofen	10 mg/mL
Beclomethasone	500 ng/mL	Mucin	1 mg/mL
Nasal Spray	5% (v/v)	Oseltamivir	500 ng/mL
Oxymetazoline	0.05 mg/mL	Phenylephrine	1 mg/mL
Zanamivir	1 mg/mL	Tobramycin	500 ng/mL

4. **Hook Effect**
Samples containing as high as 1 x 10⁷ CFU/mL of group A Streptococcus can still test positive. The tests do not show a hook or prozone effect up to the maximum observed physiological concentration (1 x 10⁷ CFU/mL).
5. **Reproducibility**
The reproducibility study is performed by two Production Technicians in production technology lab, two QC Technicians in QC lab and one R&D Technician in R&D lab respectively for four concentration levels with each concentration of sample repeated 5 times on each lot totally three lots of test devices/strips used, and tested once a day for five days. The intra-assay agreement was 100%. The inter-site agreement was 100%.

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