wellion vet

(AT) (DE) WellionVet® GLUCO CALEA Teststreifen

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(EN) WellionVet® GLUCO CALEA Test strips

IMPORTANT: PLEASE READ THIS INFORMATION AND YOUR WELLIONVET GLUCO CALEA USER GUIDE BEFORE USING THE GLUCO CALEA TEST STRIPS TO TEST YOUR BLOOD SUGAR.

INTENDED USE
The WellionVet GLUCO CALEA Blood Glucose Test Strips are to be used with the WellionVet GLUCO CALEA Blood Glucose Meters to quantitatively measure glucose in capillary whole blood.

Warnings and Precautions
[1] Do not use test strips after expiration date. The expiration date is printed on the carton and on the label of the test strip bottle. **[2]** Do not use strips that are wet or damaged. **[3]** Do not re-use the strips. **[4]** Use the WellionVet GLUCO CALEA Blood Glucose Monitoring System within temperature range of 10°C to 40°C (50°F to 104°F).

Storage and Handling
[1] Always close the vial cap tightly after removing a test strip to avoid moisture and direct sunlight. **[2]** Store the test strips at temperatures between 4°C to 30°C (39°F to 86°F). **[3]** Store out of direct sunlight; do not refrigerate or freeze. **[4]** Unopened test strips are stable until the expiration date printed on the bottle when stored properly. **[5]** The month printed refers to the END of that month. **[6]** Use within 3 months after first opening. **[7]** Do not handle the test strips with wet or dirty hands.

Test Procedure
Refer to „Testing Your Blood Glucose" section in the User's Guide of the WellionVet GLUCO CALEA Blood Glucose Meter.

Problem Solving
[1] Confirm the Test Strips are within the expiration date. **[2]** Make sure the blood completely fills the reaction zone. An Error will be displayed on the WellionVet GLUCO CALEA meter if there is insufficient blood. DO NOT add a second drop of blood; discard the test strip from the meter and retest using a new Test Strip.

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© MED TRUST, WellionVet and GLUCO CALEA are registered trademarks.