wellion®

WAVE plus

DE

HANDBUCH

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EINLEITUNG **3**

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TABLE OF CONTENTS

INTRODUCTION **26**

ABOUT BLOOD PRESSURE **27**
1. What is blood pressure? 27
2. What is hypertension and how is it controlled? 27
3. Why measure blood pressure at home? 27
4. WHO blood pressure classification 27

NOTES ON SAFETY **28**
Explanation of the icons 28
5. Blood pressure variations 28

PRECAUTIONS BEFORE USE **30**

PARTS IDENTIFICATION **30**

SYMBOLS ON DISPLAY **31**

INSERT OR REPLACE BATTERIES **31**

TIME AND SYSTEM SETUP **33**

UNIT CONVERSION mmHg/kPa DISPLAY **33**

WHO BLOOD PRESSURE CLASSIFICATION DISPLAY **34**

ATTACHING THE ARM CUFF **34**
How to take proper measurements 35

HOW TO MEASURE BLOOD PRESSURE **36**
Read memory 36
Delete Memory 36

CARE AND MAINTENANCE **37**

EN

TROUBLESHOOTING	**38**
SPECIFICATIONS	**39**
STATEMENT	**39**
FEATURE OF THE PRODUCT	**40**
REQUESTS FROM MANUFACTURER	**40**
REGARDING EMC STATEMENT	**41**
MANUFACTURER´S WARRANTY	**41**
CALIBRATION METHOD	**48**

INTRODUCTION

The Wellion WAVE plus Blood Pressure Monitor uses the oscillometric method of blood pressure measurement.

The Measurement Automatic Electronic Blood Pressure Monitor is intended for use by medical professionals or at home to monitor and display diastolic, systolic blood pressure and pulse rate, with an arm cuff around the left upper arm according to the instruction in the chapter "Attaching the arm cuff". The expected life of the product is 5 years.

The product complies with the electromagnetic compatibility requirements of IEC 60601-1-2 and safety standards of IEC 60601-1 and performance of IEC 80601-2-30 as specified in Regulation (EU) 2017/745.

ABOUT BLOOD PRESSURE

1. What is blood pressure?

Blood pressure is the force exerted by blood against the walls of the arteries. Systolic pressure occurs when the heart contracts. Diastolic pressure occurs when the heart expands.

Blood pressure is measured in millimetres of mercury (mmHg). One's natural blood pressure is represented by the fundamental pressure, which is measured first thing in the morning while one is still at rest and before eating.

2. What is hypertension and how is it controlled?

Hypertension is an abnormally high arterial blood pressure. If this is left unattended it can cause many health problems including stroke and heart attack.

Hypertension can be controlled by altering lifestyle, avoiding stress and with medication under a doctor's supervision.

To prevent hypertension or to keep it under control:


- Do not smoke
- Exercise regularly
- Reduce salt and fat intake
- Have regular physical checkups
- Maintain proper weight

3. Why measure blood pressure at home?

Blood pressure measured at a clinic or doctor's office may cause apprehension and produce an elevated reading, 25 to 30 mmHg higher than the value measured at home. Home measurement reduces the effects of outside influences on blood pressure readings, supplements the doctor's readings and provides a more accurate, complete blood pressure history.

4. WHO blood pressure classification


Standards for assessment of high blood pressure, without regard to age, have been established by the World Health Organization (WHO), as shown.
(Graph taken from Journal of Hypertension 1999, Vol 17 No.2)

EN

5. Blood pressure variations

An individual's blood pressure varies on daily and seasonal basis. It may vary by 30 to 50 mmHg due to various conditions during the day.
In hypertensive individuals, variations are even more pronounced.
Normally, the blood pressure rises during work or while playing and falls to its lowest levels during sleep. So, do not be overly concerned by the results of one measurement.
Take measurements at the same time every day using the procedure described in this manual, and know your normal blood pressure. Many readings give a more comprehensive blood pressure history.
Be sure to note date and time when recording your blood pressure. Consult your doctor to interpret your blood pressure data.

NOTES ON SAFETY

The warning signs and sample icons shown here are listed for your safety and the correct use of the unit, so as to prevent injuries or damages to the device.

The icons and meanings are as follows:

Explanation of the icons

	The icon indicates prohibitions (what you should not do). Specific actions prohibited are indicated either by text or a picture within or near the circle. The left icon refers to "general prohibition".
	The icon indicates something that is compulsory (what must always be observed). Compulsory actions are indicated either by text or pictures within or near the circle. The left icon refers to "general compulsion".
	Do not disassemble

	Type BF applied part Electrically connected to patient but not directly to heart
	Caution
	Please refer to the instructions for use
	Consult instructions for use
	Contact your local authorities to determine the proper method for disposal of potentially bio hazardous parts and accessories.

	Indicates a medical device that needs to be protected from moisture.
MD	Medical device
IP20	Ingress protection degree as defined by IEC 60529 IP20: Product is touchproof and will be resistant to dust or objects that are over 12mm in size (e.g. persons fingers), but is not waterproof.

	MR Unsafe An item marked MR unsafe is known to pose hazards in all MR environments.
	Importer


Only a healthcare professional is qualified to interpret blood pressure measurements. Self-measurement means control, not diagnosis or treatment. Unusual values must always be dis-cussed with your doctor. The device should not be used to decide on first aid measures nor for continuous monitoring of blood pressure. This device is not suitable for measurement in any case of medical emergency or during operations. Although portable, it should not be operated in transit (e.g. patient transport). It is intended for use at home or in any stationary location. Please press the "ON/OFF" button to stop the device from working when you feel uncomfortable with the arm, or if the air is inflating abnormally without stop.	Caution	
The device should not be operated by children under 18 or people with restricted physical, sensory, or mental acumen so to avoid hazardous situations.		
Do not use the device for purposes other than measuring blood pressure, as it may cause an accident or trouble.		
During measurement, make sure that no mobile phone or any other electrical devices emitting electromagnetic fields is in the surrounding area. This may result in incorrect operation of the device and/or cause an inaccurate reading.		
The device is prohibited from being used during movement.		
Do not use the equipment outdoors or in wet areas		
Stop using the equipment immediately, if it is in contact with water.		
Do not disassemble, repair, or remodel the main unit or the cuff of the blood pressure monitor. This will lead to malfunctioning of the device		

EN

PRECAUTIONS BEFORE USE


1. If you are taking medication, consult with your doctor to determine the most appropriate time to measure your blood pressure. NEVER change a prescribed medication without first consulting with your doctor.
2. For people with irregular or unstable peripheral circulation problems due to diabetes, liver disease, hardening of the arteries, etc., there may be fluctuation in blood pressure values measured at the upper arm versus at the wrist.
3. Measurements may be impaired if this device is used near televisions, microwave ovens, X-ray, mobile phone equipment or other devices with strong electrical fields.
 To prevent such interference, use the monitor at a sufficient distance from such devices or turn them off.
4. Wash your hands before use.
5. Do not measure on the same arm which uses simultaneously another medical monitoring ME Equipment, otherwise it could cause loss of function.
6. Consult your doctor if unexpected readings are obtained, also please refer to "TROUBLESHOOTING" of the manual.
7. As people generally are more relaxed at home than at the doctor´s office, blood pressure may be lower at home.
8. Cuff pressure range 0-299 mmHg

PARTS IDENTIFICATION


EN


SYMBOLS ON DISPLAY

Accessory


Manual

Adapter

INSERT OR REPLACE BATTERIES

1. Remove the battery cover.
2. Insert new batteries into the battery compartment as shown, taking care that the polarities (+) and (-) are correct.
3. Close the battery cover. Use only LR6, AA batteries.

1.5V batteries
(LR6 or AA)


Disposal of empty battery to the authorized collecting party subject to the regulation of each individual territory.


EN

⚠ Caution:
- Insert the batteries as shown in the battery compartment. If not, the device will not work.
- When „LOW BATTERY“ symbol blinks in the display, replace all batteries with new ones. Do not mix old and new batteries. It may shorten the battery life, or cause the device to malfunction.
 „LOW BATTERY“ symbol does not appear when the batteries run out.
- Please ensure to match the correct polarity when replacing batteries.
- Batteries, which have a fluid on the surface or have been modified, must not be inserted into the products.
- Battery short circuit must be prevented.
- Battery lifespan varies with the ambient temperature and may be shorter at low temperatures.
- The batteries may leak and cause a malfunction.
- Use only batteries corresponding to given specifications. The batteries provided with the device are for testing monitor performance and may have a shorter life.
- Used batteries may leak and damage the main unit. Please observe the following points:
 - Remove batteries if the device is not likely to be used for some time (approximately three months or more).

TIME AND SYSTEM SETUP


1. Press the „SET“ button to enter the time setting menu.
2. To enter the time setting mode when switched off, press and hold the “SET” button until the year flashes on the LCD.
3. Press the “MEM” button to adjust the year, then press the “SET” button again to save your settings and enter the month setting mode.
4. Press the “MEM” button to adjust the month. Repeat the same steps to adjust day/hour/minute until setting is completed.

UNIT CONVERSION mmHg/kPa DISPLAY

Blood pressure readings can be displayed in either one of the following units: mmHg or kPa (mmHg being the preset unit).
To enter the unit settings menu press the „ON/OFF“ button for 10 sec. Choose the unit of your liking by pressing the „MEM“ button and then exit the menu by pressing the „ON/OFF“ button again. The current unit is displayed on the LCD screen.
Please note that all memorized results change their unit as well.

WHO BLOOD PRESSURE CLASSIFICATION DISPLAY

Diastolic blood pressure Reference material:
Journal of hypertension 1999. vol 17 No.2


- **Grade 3 hypertension** (severe)
- **Grade 2 hypertension** (moderate)
- **Grade 1 hypertension** (mild)
- **High-normal**
- **Normal**
- **Optimal**

ATTACHING THE ARM CUFF


1. Wrap the arm cuff around the upper arm, about 2-3 cm above the elbow, as shown. Place the cuff directly on the skin, as clothing may cause a faint pulse and result in a measurement error.
2. Constriction of the upper arm, caused by rolling up a shirtsleeve, may prevent accurate readings.
3. Secure the arm cuff with the Velcro Strip in such a way that it lies comfortably and is not too tight. Lay the arm on the table (palm upwards) so that the arm cuff is at the same height as the heart. Make sure that the tube is not kinked.
4. Choose the appropriate cuff size dependent on your arm circumference (refer to "Specifications").

Avoid constriction of the upper arm

How to take proper measurements


For best accuracy in blood pressure measurement:

- Patient in seated position, including
 1) comfortably seated
 2) legs uncrossed
 3) feet flat on the floor.
 4) back and arm supported
 5) middle of the CUFF at the same level as the heart
- Remain still and keep quiet during measurement.
- Relax as much as possible and do not talk during the measurement process.
- Measure your blood pressure at about the same time every day.
- Do not measure right after physical exercise or after a bath. Take a rest for twenty or thirty minutes before taking the measurement.
- Please note that the following conditions can influence the measurement result:
 - eating a meal, consuming drinks containing alcohol or caffeine, exercising or bathing up to one hour before the measurement is taken
 - talking, nervousness, mood swings, drastic changes in ambient temperature during the measurement
 - measurement in moving vehicles
 - continuous measurement over a longer period of time

EN

HOW TO MEASURE BLOOD PRESSURE

1. Set up the arm cuff to your upper arm as described in the previous section of "ATTACHING THE ARM CUFF".
2. Press the "ON/OFF" button, all icons automatically appear for two seconds on the display, then the device switches to measurement and displays "00" or last measurement record.

3. When the measurement starts, the cuff in the strap will automatically inflate. This is indicated by the "Inflation" symbol "♥". When the measurement is completed, the results will be displayed on the LCD display.


Read memory


Press the "MEM" button while the meter is off to retrieve saved data. The meter shows the average value of the latest 3 measurements on the display, indicated by the „ AVG" symbol in the right upper corner of the display. In order to retrieve your individual blood pressure values, press the „MEM" button again. To go through the memory, press „MEM" for navigating up or „SET" for navigating down. It is also possible to review your measurements directly after measurement by pressing the „MEM" button.


Delete Memory

Enter the "READ MEMORY" mode, press and hold the "MEM" button for five seconds until the LCD displays "no". This indicates that all stored results have been deleted successful.

CARE AND MAINTENANCE

To keep your digital blood pressure monitor in the best condition and to protect the unit from damage, follow the directions listed below:

- Keep the monitor in the storage case when not in use.
- Do not fold the arm cuff too tightly.
- Clean the monitor and cuff with a soft dry cloth.
- Do not use any abrasive or volatile cleaners.
- Cuff shall be detached from the device softly and slowly.
- Cleaning of the cuffs is recommended after each use.
 - Use a soft cloth lightly moistened with 75% ethanol to wipe the entire inner side of the cuff which is in contact with the skin. Wiping shall be done for a total of 3 times, each time using a fresh cloth.
 - Air-dry the cuff.

Caution:
- Do not submerge the device or any of the components in water.
- Do not subject the monitor to extreme hot or cold temperatures, humidity or direct sunlight.
- Store the device and the components in a clean, safe location.
- Do not drop this monitor or subject it to strong impact.
- Remove the batteries if the unit will not be used for three months or longer.
- Always replace the whole lot of batteries at the same time.

Information:
- We recommend calibration checks every two years for devices used professionally, for example in pharmacies, doctors' offices, or clinics. Please also observe the national regulations specified by law. Calibration checks may only be performed by the competent authorities or authorized maintenance services.
- For lay use, the device complies with the permitted maximum measurement deviation throughout its entire service life.

See the Calibration Method for more details.

EN

TROUBLESHOOTING

If you have trouble in using the product please check the following points first:

Error display	Possible cause	How to correct
Device doesn't power on when ON/OFF button is pressed or the battery icon flashes	No batteries installed.	Insert batteries
	Batteries completely drained.	Replace batteries
	Polarity of batteries is not matched correctly.	Make sure to match the polarity of the batteries correctly.
E0: Device can`t detect the signal pulse	Cuff is worn over clothing.	The cuff should preferably be worn on the bare upper arm or only over thin clothing without rolling up your sleeves.
	The cuff may be worn in wrong direction and the tube interface doesn`t correspond to the arterial position.	
	If the E0 error code is not caused by the above mentioned reasons the device is defect.	
E1: cuff cannot be inflated properly	Check your arm cuff for any air leakages.	Replace arm cuff with a new one
E3: inflating pressure is too high	Pressure value of more than 299mmHg	Start a re-measurement or return the device for recalibration.
E2E4: shaking while measurement	Hand or body shaking while measurement	Keep correct position and gesture to measure again
Battery icon displayed	Batteries have low power	Replace batteries and repeat the measurement
The systolic or diastolic pressure value is too high	The arm cuff was held lower than your heart level	Keep correct position and gesture to measure again
	The arm cuff was not attached properly	
	You moved your body or spoke during measurement	
The systolic or diastolic pressure value is too low	The arm cuff was held higher than your heart level	
	You moved your body or spoke during measurement	

SPECIFICATIONS

Measuring Method:	Oscillometric Measurement	**Indication:**	Digital LCD display
Measuring Range:	Pressure: 30 - 280 mmHg Pulse: 40 - 199 beats/minute	**Accuracy:**	Static Pressure: ±3mmHg Pulse: ± 5%
Memory:	Stores up to 90 blood pressure readings	**Power supply:**	4x 1,5V batteries (LR6 or AA) approximately 200 measurements (using new alkaline batteries)
Operating condition:	+ 5°C - + 40°C; 15% RH - 93% RH; Atmospheric pressure: 70kPa - 106kPa	**Storage condition:**	- 20°C - + 55°C; 0% RH - 93% RH; Atmospheric pressure: 50kPa - 106kPa
Dimensions:	Approx: 140(W) X 110(H) X 48D) mm	**Weight:**	Approx: 410 g, excl. batteries
Classification:	Type BF	**Upper arm circumference:**	22 - 42 cm
AC adapter:	Input: 100 – 240V 50/60Hz, 0,2A Output: 6V, 500 mA		

1. Type of protection against electric shock: INTERNALLY POWERED EQUIPMENT.
2. Degree or protection against electric shock: TYPE BF APPLIED PART.
3. Mode of operation: CONTINUOUS OPERATION.
4. Equipment not suitable for category AP & APG equipment use in presence.

STATEMENT

The system might not meet its performance specifications if stored or used outside the temperature and humidity as mentioned below:

- Operating conditions: +5°C - +40°C 15%RH - 93%RH; 70kPa - 106kPa
- Storage conditions: -20°C - +55°C 0%RH - 93%RH

EN

FEATURE OF THE PRODUCT

1. Stores up to 90 blood pressure readings.
2. Large and clear LCD display.
3. WHO blood pressure classification display.
4. Easy to use: Press a single button to start automatic measurement. Both measurement values and measurement time are recorded.
5. Automatically turns off (within 1 minute) to save power.

REQUESTS FROM MANUFACTURER

- Make sure that the connection tube is not kinked before starting the measurement to avoid any injury to patient.
- Do not take blood pressure measurements more often than three times in a row. In between measurements there should be a resting interval of at least 5 minutes in order to prevent blood flow interference.
- DO NOT take blood pressure measurements more often than necessary, do not exceed a number of 6 measurements per day.
- Do not apply the cuff over a wound as this can cause further injury.
- Do not measure on the arm which is on the side of a mastectomy, otherwise it could cause injury.
- Observe the air pressure value from the LCD display.
- When measuring, the device cannot exceed 280 mmHg, otherwise please press the "ON/OFF" button to stop.
- Do not use force to bend the arm cuff or the air tube.
- DO NOT drop or subject the device to severe shocks or vibrations.
- Please use only accessories and detachable parts specified/ authorised by the manufacturer. Otherwise, it may cause damage to the unit or pose danger to the user/patients.
- For service information, parts list etc., please contact the dealer.

- The PATIENT is an intended OPERATOR.
- The device must only be serviced and maintained when switched off.
- The device may be maintained by the user itself as described in the according section of this manual.
- Stop using the equipment immediately after contact with water.

REGARDING EMC STATEMENT

1. A statement to the effect that this equipment needs to be installed and put into service in accordance with the information provided in the ACCOMPANYING DOCUMENTS; and
2. a statement to the effect that wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect this equipment and should be kept distanced away from the equipment. The distance is calculated by the MANUFACTURER from the 800 MHZ to 2,5 GHz column of Table 5 or Table 6 of EN 60601-1-2:2015, as appropriate.

MANUFACTURER´S WARRANTY

The producer warrants Wellion WAVE plus blood pressure meter against defects in materials and workmanship for a period of two years from the date of purchase. The warranty is lost if the instrument is misused, poorly maintained or is tampered with. Liability under this warranty is limited to the repair of defective parts or - at the discretion of the producer - to the replacement of the instrument. The right to rescind the purchase agreement exists only if the replacement is also found to be defective. Claims other than these will not be dealt with. The warranty is not valid if the damage results from misuse, maltreatment, tampering, human error and the use of extreme force. This warranty is valid only if the date, the stamp and the signature of the dealer are recorded on the warranty card on the date of purchase. The warranty period is not to be extended by any claim made under this warranty.

EN

Appendix 1 Guidance and Manufacturer Declaration Tables

Guidance and manufacturer´s declaration- electromagnetic emission		
The Model PG-800B11 is intended for use in the electromagnetic environment specified below. The customer or the user of the Model PG-800B11 should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Model PG-800B11 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Model PG-800B11 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	A	
Voltage fluctuations /flicker emissions IEC 61000-3-3	Complied	

Guidance and manufacturer´s declaration – electromagnetic immunity			
The Model PG-800B11 is intended for use in the electromagnetic environment specified below. The customer or the user of the Model PG-800B11 should assure that it is used in such an environment.			
Immunity test	**IEC 60601 Test level**	**Compliance level**	**Electromagnetic environment -guidance**
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	± 8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrostatic transient/ burs IEC 61000-4-4	± 2 kV for power supply lines 100 kHz repetition frequency ±1 kV for input/output lines	± 2 kV for power supply lines 100 kHz repetition frequency ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±0.5 kV, ±1kV differential mode line-line	±0.5 kV, ±1kV differential mode line-line	Mains power quality should be that of a typical commercial or hospital environment
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT (100% dip in UT) for 0.5 cycle at 0°, 45°, 90°, 135°,180°, 225°, 270° and 315° 0% UT (100% dip in UT) for 1 cycle at 0° 70% UT (30%dip in UT) for 25/30 cycles at 0° 0% UT (100% dip in UT) fpr 250/300 cycle at 0°	0% UT (100% dip in UT) for 0.5 cycle at 0°, 45°, 90°, 135°,180°, 225°, 270° and 315° 0% UT (100% dip in UT) for 1 cycle at 0° 70% UT (30%dip in UT) for 25/30 cycles at 0° 0% UT (100% dip in UT) fpr 250/300 cycle at 0°	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Model PG-800B11 product name requires continued operation during power mains interruptions, it is recommended that the Model PG-800B11 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m,50/60Hz	30 A/m,50/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: UT is the a. c. mains voltage prior to application of the test level.			

EN

Guidance and manufacturer's declaration-electromagnetic immunity			
The Model PG-800B11 is intended for use in the electromagnetic environment specified below. The customer or the user of the Model PG-800B11 should assure that it is used in such an electromagnetic environment.			
Immunity test	**IEC 60601 test level**	**Compliance level**	**Electromagnetic environment - guidance**
Conducted RF IEC 6100-4-6	3 Vrms 150 kHz to 80 MHz 6 Vrms 150 kHz to 80 MHz outside ISM bandsa	6 V	Portable and mobile RF communications equipment should be used no closer to any part of the Models PG-800B11, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. **Recommended separation distance** $d = \left[\frac{3.5}{V_1}\right]\sqrt{P}$ $d = \left[\frac{12}{V_2}\right]\sqrt{P}$ 80MHz to 800MHz $d = \left[\frac{12}{E_1}\right]\sqrt{P}$ 800MHz to 2.7GHz $d = \left[\frac{23}{E_1}\right]\sqrt{P}$

Radiated RF IEC 6100-4-3	10 V/m 80 MHZ to 2.7 GHz	10 V/m	Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres(m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range b Interference may occur in the vicinity of equipment marked with the following symbol:

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.

b The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,7 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in these frequency ranges.

c Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Model PG-800B11 is used exceeds the applicable RF compliance level above, the Model PG-800B11 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Model PG-800B11.

d Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

EN

Recommended separation distances between portable and mobile RF communications equipment and the Model PG-800B11				
The Model PG-800B11 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Model PG-800B11 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Model PG-800B11 as recommended below, according to the maximum output power of the communications equipment.				
Rated maximum output of transmitter	Separaton distance according to frequency of transmitter (m)			
	150 kHz to 80 MHz ISM bands $d = \left[\frac{3.5}{V_1}\right]\sqrt{P}$	150 kHz to 80 MHz ISM bands $d = \left[\frac{12}{V_2}\right]\sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{12}{E_1}\right]\sqrt{P}$	800 MHz to 2.7 GHz $d = \left[\frac{23}{E_1}\right]\sqrt{P}$
0.01	0.12	0.12	0.12	0.23
0.1	0.38	0.38	0.38	0.73
1	1.2	1.2	1.2	2.3
10	3.8	3.8	3.8	7.27
100	12	12	12	23

For transmitters rated at a maximum output power not listed above the recommended separation distance *d* in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0, 15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.

NOTE 3: An additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,7 GHz to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas.

NOTE 4: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

EN

CALIBRATION METHOD

1. Press and hold the „ON/OFF, MEM“ button at the same time, load the battery, enter the static air pressure mode after the LCD screen is fully displayed, and then release the button.
2. Press ON/OFF to close the internal air valve.
3. Connect the external standard barometric interface and the digital barometer interface to the cuff interface.
4. External input 50mmHg and 200mmHg standard static air pressure, and observe the air pressure value displayed at the position of the LCD systolic pressure (SYS) and the value of the digital pressure gauge should be in the range of +/- 3mmHg.

⚠ Caution:
1. ME devices can be used in exposed environments, including electromagnetic interference environment to ensure basic safety and basic performance unchanged.
2. In the event of any serious event related to this product, such as serious adverse event, significant alteration of the product resulting in change of intended use, etc., it will be reported to the manufacturer and the competent authorities of the user and/or the member states where the patient is located.

Notes:
Essential performance: Limits of the error to the manometer, ±3mmHg. Reproducibility of the blood pressure determination, ±3mmHg.
Clinical benefits: Accurate measurement of SBP and DBP, clinical performance meets the requirements of ISO 81060-2:2018.

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Manufacturer: Shenzhen Pango Medical Electronics Co.,Ltd
Tel: 86-755-33825988, Fax: 86-755-33825989
Main Site: Building 2, No. 25 Fenghuang Road, Industrial Zone, Xikeng First Village, Henggang Street, Longgang District, Shenzhen, 518115 Guangdong, P. R. China.
Additional site1:2-4 Floor ,No.5 Shanzhuang Rd., Xikeng Village, Henggang Town, Longgang District, Shenzhen, 518115 Guangdong, P. R. China.

MED TRUST
THE MEDICAL SERVICES COMPANY

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MD **MEDICAL DEVICE**

WWAVE032PLUS 20260113
Model No.: PG-800B11
Software Version 2.5

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© MED TRUST, Wellion and WAVE are registered trademarks.
Lotus NL B.V.